Artritest

directo

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1 x 0,5 ml Controle Negativo

REFERÊNCIA

- Singer, J.M. and Plotz, C. - Am. J. Clin. Path. 28/611 (1957).
- Singer, J.M. et al. - Am. J. Clin. Path. 72/4:597 (1979).
- Brooks, G.W. and Cobb, S. - Arthritis and Rheumatism 6/3:198 (1963).
- Plotz, C.M. and Singer, J.M. - Annual Meeting of the American Rheumatism Association, Chicago (1956).
- Oreskes, I.; Singer, J.M. and Plotz, C.M. - J. Immunology 90:107 (1963).
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 4th ed., 2001.

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Artritest

directo

Direct agglutination slide test for the diagnosis of
rheumatoid arthritis

SUMMARY

Rheumatoid arthritis is a chronic syndrome of unknown etiology, characterized by an unspecific and generally symmetrical inflammation of the joint, which sometimes evolves towards the destruction of the joint and periarticular structures.

The affected joint undergoes a thickening of the synovial membrane with the formation of pleats and lymphocytes proliferation. These cells, which form lymphoid follicles, are responsible for the Rheumatoid Factor (RF) and other immunoglobulins. RF are generally, IgM anti-IgG class. They react with IgG Fc fractions. RF of IgG class have also been found, but in a lesser proportion of cases.

RF class are present in 85-90 % of adults suffering from rheumatoid arthritis. However, they are not specific of rheumatoid arthritis. They have also been found in 3-5% healthy individuals.

For this reason, when diseases other than rheumatoid arthritis are suspected, and there is a positive agglutination reaction, it is advisable to perform globulins precipitation with ammonium sulfate. In this manner, RF is thoroughly precipitated and are kept in solution any agglutinin, responsible for a false positive reaction, as well as any possible inhibitor.

PRINCIPLE

The IgM class rheumatoid factor is detected in presence of gamma-immunoglobulin or Cohn's fraction II (in this case the antigen) adsorbed onto an inert latex polystyrene support. This antigen binds to RF (anti-IgG), producing latex polystyrene particles agglutination, macroscopically visible.

PROVIDED REAGENTS

A. Reagent A: suspension of 0.20 microns latex polystyrene particles, which have thermoaggregated molecules of gamma-immunoglobulin or Cohn's fraction II adsorbed.

Positive Control: dilution of serum proteins containing rheumatoid factor. Equivalent to 2 (+).

Negative Control: dilution of non-reactive serum proteins.

INSTRUCTIONS FOR USE

Reagent A: shake well before use, emptying the dropper pipette first.

Positive and Negative Controls: ready to use.

WARNINGS

Reagents are for "in vitro" diagnostic use.

The controls have been tested and found non-reactive to HIV, HCV and HBV. However, they should be handled as infectious material.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: are stable at 2-10°C until the expiration date shown on the box. Do not freeze.

INSTABILITY OR DETERIORATION OF REAGENTS

Self agglutination of Reagent A is a sign of its deterioration; in that case, discard.

SAMPLE

Serum

a) Collection: obtain serum in the usual way. Do not use heat-inactivated samples.

b) Additives: not required.

c) Known interfering substances: hemolysis and certain proteins (other than the rheumatoid factor) may cause erroneous results.

See Young, D.S. in References for effect of drugs on the present method.

d) Stability and storage instructions: serum should be preferably fresh. If it cannot be processed immediately, it can be keep at 2-10°C no longer than a week from its collection.

REQUIRED MATERIAL

1- Provided

- plastic or glass slides with black background
- 25 ul dropper cap

2- Non-provided

- Wood or glass rods
- Stopwatch
- Pipettes and micropipettes for measuring the stated volumes
- Light source

PROCEDURE

Bring the reagents and samples to room temperature before use.

I- QUALITATIVE TECHNIQUE

In one of the circles of the provided slide, place:

Sample or Control	25 ul
Reagent A	1 drop (25 ul)

Mix with a wood or glass rod to obtain a uniform suspension all over the circle surface and immediately start the stopwatch. Gently move the slide and observe the result macroscopically under a light source, during 2 minutes.

II- TITRATION

Positive sera can be titrated performing serial dilutions, in 6 Kahn tubes.

1) Place 0.4 ml saline solution in tube N° 1 and 0.2 ml saline solution in the other 5 ones.

2) Add 0.1 ml (100 ul) serum to tube N° 1 and mix. Transfer 0.2 ml of this dilution to tube N° 2 and mix, continuing with the dilutions in this manner up to the last tube. The obtained dilutions are 1:5; 1:10; 1:20; 1:40; 1:80; 1:160.

3) Test each dilution following the Qualitative Technique.

without evident symptoms of rheumatic disease, 72 non-reactive samples were obtained. Therefore, it was concluded that 4% of healthy individuals are RF positive.

WIENER LAB. PROVIDES

Kits for 150 tests (Cat. N° 1103153) containing:

1 x 4 ml Reagent A

1 x 0.5 ml Positive Control

1 x 0.5 ml Negative Control

REFERENCES

- Singer, J.M. and Plotz, C.M. - Am. J. Clin. Path. 28/611 (1957).
- Singer, J. M. et al. - Am. J. Clin. Path. 72/4:597 (1979).
- Brooks, G. W. and Cobb, S. - Arthritis and Rheumatism 6/3:198 (1963).
- Plotz, C M. and Singer, J.M. - Annual Meeting of the American Rheumatism Association, Chicago (1956).
- Oreskes, I.; Singer, J.M. and Plotz, C.M. - J. Immunology 90:107 (1963).
- Young, D. S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 4th ed., 2001.

INTERPRETATION OF RESULTS

Negative: homogeneous suspension.

Positive: agglutination that appears within 2 minutes. It is qualified from 1 to 4 (+).

Titer: is the inverse of the highest dilution showing a macroscopically visible agglutination.

The approximate rheumatoid factor concentration of the sample can be calculated with the following formula:

$RF \text{ (IU/ml)} = \text{Titer} \times \text{Reaction Sensitivity (1 IU/ml)}$

E.g.: the sample has a 1:40 titer. Its RF concentration is $40 \times 1 = 40 \text{ IU/ml}$.

QUALITY CONTROL METHOD

Simultaneously process the provided Controls or tested reactive serum, using 25 ul of the corresponding Control and 25 ul of Reagent A following the qualitative technique.

REFERENCE VALUES

They have not been clearly established. Values above 20-30 IU/ml are considered pathological.

PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE.

Any variation in the Sample/Reagent proportion may produce erroneous results.

PERFORMANCE

a) Analytic sensitivity: calibrated against the CDC Standard, a sensitivity of 1 IU/ml is obtained.

In a study performed over 176 samples from a hospitalized population, the results were compared with a quantitative turbidimetric method, obtaining a sensitivity of 100% and a specificity of 99.3%.

In another study performed on a population of 33 samples from patients with a medical record of rheumatoid arthritis, 29 reactive samples were obtained. It was concluded that 88% of patients with rheumatoid arthritis are RF positive.

b) Specificity: in a study performed over 75 samples from 50 men and 26 women aged between 18 and 47 years old,

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Wiener Laboratorios S.A.I.C.