Chagatest

ELISA recombinante v.3.0

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Chagatest

ELISA recombinante v.3.0

3rd Generation method for the detection of antibodies
to *Trypanosoma cruzi*

SUMMARY

Chagas' disease is a parasitic infection caused by the *Trypanosoma cruzi*.

The laboratory diagnosis depends on the stage of the disease. During the acute phase, the diagnosis is directly done by detecting the parasite in blood or through immunological methods that detect IgM. During the chronic phase, immunological methods can be used, such as the Machado and Guerreiro complement fixation reaction or latex agglutination, flocculation, hemagglutination, direct agglutination, immunofluorescence and lately ELISA.

PRINCIPLE

In this qualitative technique for the detection of antibodies anti-*T. cruzi*, the sample is diluted in the support in which the recombinant antigen (1, 2, 13, 30, 36 and SAPA) is immobilized (3rd generation method). These antigens are obtained by DNA recombinant techniques starting from specific proteins from the epimastigote and trypomastigote stages of the *T. cruzi* corresponding to highly conserved zones among different strains. The technology used allows us to ensure an antigenic mixture of known and permanent composition batch to batch, giving reproducible, specific and highly sensitive results. If the sample contains specific antibodies, these will form a complex with the antigens and will remain bound to the support. The unbound fraction is eliminated by washing, after which anti-human immunoglobulin antibodies conjugated to peroxidase are added. If a reaction was produced in the first step of the process, the conjugate will be bound. After a new wash the enzymatic substrate is added. If bound conjugate is present, a light blue color will be developed. The reaction is stopped by adding sulfuric acid, and the color change to yellow.

PROVIDED REAGENTS

Coated microtitration plate: microtitration plate with removable strips containing immobilized *Trypanosoma cruzi* recombinant antigens.

Conjugate: goat anti-human immunoglobulins conjugated to peroxidase.

Substrate A: 60 mmol/l hydrogen peroxide in 50 mmol/l citrate buffer, pH 3.2.

Substrate B: 0.01 mol/l tetramethylbenzidine (TMB) in 0.1 N hydrochloric acid.

Stopper: 2 N sulfuric acid.

Concentrated Wash Buffer: 1.4 mol/l sodium chloride in 100 mmol/l phosphate buffer and 0.1 g/l non-ionic surfactant.

Sample Diluent: bovine albumin in phosphate buffered saline (PBS), pH 7.2.

Positive Control: dilution of inactivated serum containing antibodies against *Trypanosoma cruzi*.

Negative Control: dilution of inactivated non-reactive serum.

INSTRUCTIONS FOR USE

Wash Buffer: to use dilute 1+ 4 with distilled water (1 part Wash Buffer + 4 parts distilled water).

Constituents of this reagent may precipitate at low temperature. In such case, place in water bath at 37°C for a few minutes and mix by inversion.

Coated microtitration plate, Conjugate, Substrate A, Substrate B and Stopper: ready to use.

Sample Diluent: ready to use. The Sample Diluent color may vary from batch to batch without affecting the performance.

Negative and Positive Controls: ready to use.

WARNINGS

- All samples should be handled as capable of transmitting infection. Controls are inactivated, however, they should be handled as infectious material.

- The control sera have been examined for Hepatitis B Surface Antigen (HBsAg), Human Immunodeficiency Virus (HIV) and Hepatitis C virus and have been found non-reactive.

- All the materials used to perform the test must be destroyed to ensure the inactivation of pathogenic agents. The recommended method is autoclaving one hour at 121°C. Liquid waste can be disinfected with sodium hypochlorite (5% final concentration) for at least 60 minutes.

- Do not exchange reagents from different kits and batches.

- Do not use foreign reagents.

- The microtitration plates should be incubated in incubator. Do not use water bath. Do not open the incubator during this process.

- Avoid the contact of strips with hypochlorite fumes, from biohazardous disposal containers or other sources, because hypochlorite affects the reaction.

- Reagents are for "in vitro" diagnostic use.

- Sulfuric acid is corrosive. H315+H320: Causes skin and eye irritation. H314: Causes severe skin burns and eye damage.

P262: Do not get in eyes, on skin, or on clothing. P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P302 + P352: IF ON SKIN: Wash with plenty of soap and water. P280: Wear protective gloves/protective clothing/eye protection/face protection.

- Avoid the spill of liquids and the formation of sprays.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents are stable in refrigerator (2-10°C) until the

expiration date shown on the box. Do not freeze.

Wash Buffer: stable for 3 months at room temperature.

Coated microtitration plate: well strips with immobilized antigen are provided in vacuum-sealed pouches with desiccant. Do not open the pouch until performing the test and it has reached room temperature, otherwise the contents can get moistened. Unused well strips should be kept in the pouch with the desiccant, sealed and stored at 2-10°C. Under these conditions, the strips can be used within 5 months, as long as the expiration date of the kit has not elapsed.

SAMPLE

Serum or plasma

a) Collection: obtain in usual way. Do not use samples inactivated by heat. See PROCEDURE LIMITATIONS.

b) Additives: not required for serum. When using plasma, any anticoagulant commonly used for transfusions can be employed.

c) Known interfering substances: hemolysis, hyperlipemia or other causes of turbidity may cause erroneous results. These samples must be clarified by centrifugation.

d) Stability and storage instructions: undiluted samples can be stored at 2-10°C for 7 days. For longer storage periods samples should be frozen at -20°C or lower. Avoid repeated freezing and thawing. Samples should not be refrozen, since it may cause erroneous results.

If samples are to be transported, they should be packed according to local regulations for biohazardous materials shipment.

REQUIRED MATERIAL (non-provided)

- Micropipettes for measuring the stated volumes.
- Timer or stopwatch.
- Incubator at 37°C.
- Spectrophotometer for microtitration plates reading (optional).

ASSAY CONDITIONS

- Primary wavelength: 450 nm
- Secondary wavelength (bichromatic): 620-650 nm
- Instrument calibration: bring the spectrophotometer to zero with Reagents Blank, processing it as a test but omitting addition of the Sample.
- Total reaction time: 90 minutes
- Reaction temperature: 37°C and room temperature
- Sample volume: 10 µl

PROCEDURE

Bring the reagents and the samples to room temperature before use. Once the test has begun, it should be completed without interruption.

Simultaneously process 2 Positive Controls (PC), 3 Negative Controls (NC) and the unknown (U). When placing the sample and/or Controls on the Sample Diluent, make sure to put them on the center of the liquid and not on the walls or bottom of the well. Rinse the pipette with the Diluent

dispensed in the well to ensure a correct homogenization. In the microtitration plate to be used, place:

	U	PC	NC
Sample Diluent	200 µl	200 µl	200 µl
Positive Control	-	10 µl	-
Negative Control	-	-	10 µl
Sample	10 µl	-	-

Once the samples have been dispensed on each strip, mix by gently tapping the sides of the microtitration plate for 10 seconds. To avoid evaporation, cover the plate and incubate 30 minutes in incubator at 37°C. Then, carefully aspirate the fluid from each well into a biohazardous container with 5% sodium hypochlorite. Then, wash 5 times with Wash Buffer, using approximately 300 µl/time/well. After each wash, the liquid should be discarded into the container with hypochlorite. Optionally use automated washer. After the last wash, completely discard the residual liquid, inverting and tapping the microtitration plate several times onto absorbent paper, gently pressing with the hand the long sides to prevent the well strips from falling. Then, add to each well:

Conjugate	1 drop	1 drop	1 drop
------------------	--------	--------	--------

If using an automated micropipette, dispense 50 µl. Mix by gently tapping the sides of the microtitration plate for 10 seconds. To avoid evaporation, cover the plate and incubate 30 minutes in incubator at 37°C. Then, carefully aspirate the fluid from each well into the container with hypochlorite, and wash as indicated above. After the last wash, completely discard the residual liquid, inverting and tapping the microtitration plate several times onto absorbent paper, gently pressing with the hand the long sides to prevent the well strips from falling. Then, add to each well:

Substrate A	1 drop	1 drop	1 drop
Substrate B	1 drop	1 drop	1 drop

If using an automated micropipette, dispense 50 µl. Mix by gently tapping the sides of the microtitration plate for 10 seconds. Incubate 30 minutes at room temperature. Then add:

Stopper	1 drop	1 drop	1 drop
----------------	--------	--------	--------

If using an automated micropipette, dispense 50 µl. Mix by gently tapping the sides of the microtitration plate for 10 seconds. Read in spectrophotometer at 450 nm or bichromatic at 450/620-650 nm or evaluate the result with the naked eye by comparing it with Positive and Negative Controls.

STABILITY OF FINAL REACTION

Reaction color is stable for 30 minutes, therefore the results should be read within this period.

RUN VALIDATION CRITERIA

Run is valid if the following two conditions are met simultaneously:

a) The readings of at least 2 out of the 3 Negative Controls corrected against the Reagents Blank, should be equal or less than 0.150 O.D.

b) The corrected mean reading of Positive Controls should be equal or higher than 0.600 O.D.

If one or both of these conditions are not met, repeat the run. Remember that the obtained readings will depend on the sensitivity of the reader in use.

INTERPRETATION OF RESULTS

a) With optical instruments:

The presence or absence of antibodies to *T. cruzi* is determined by relating the sample absorbance with the Cut-off value.

Cut-off = NC + 0.300 O.D.

where NC: average of Negative Controls readings

Indetermination Area: Cut-off $\pm$ 10%

Non-Reactive Samples: samples with absorbances lower than the indetermination area lowest limit.

Reactive Samples: samples with absorbances higher than the indetermination area highest limit.

Indeterminate Samples: samples with absorbances falling within the indetermination area. These samples must be retested.

b) Visual interpretation:

If this kind of interpretation is chosen, all samples that do not show a stronger color than the Negative Controls, should be considered non-reactive. On the contrary, a yellow sample is considered reactive.

PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE.

- The following are causes of erroneous results:
 - Improper washing of wells.
 - Cross contamination of non-reactive samples with antibodies from a reactive sample.
 - Contamination of the chromogenic solution with oxidizing agents (chlorine, etc.).
 - Contamination of the Stopper.
 - Improper storage of unused strips wells.
 - Do not use water bath for the incubation.
 - Contamination of the diluted Wash Buffer. Verify that clean bottles are used to prepare and store the diluted Wash buffer. If turbidity or precipitation are observed when prepared, it should be discarded.
- A negative result does not exclude the possibility of exposure or infection by *T. cruzi*.
- Occasionally negative absorbances may be obtained when performing bichromatic readings that do not invalidate the determination. This is because some samples yield readings below the Reagent Blank.
- Samples that may yield values above the numeric reading range are found, depending on the sensitivity of the photometer used. These results should be interpreted as reactive. For these particular cases reading at 490 nm or 490/650 nm can be performed in order to obtain numeric values.

- Do not use heat inactivated samples, since they may yield false positive results.
- Check that the washing system in use (WIENER WASHER or other) aspirates properly, and that the wash solution is dispensed evenly.

PERFORMANCE

a) Sensitivity: over a panel of 70 samples with positive xenodiagnostic and serology, the sensitivity is 100%.

Over another panel of 144 samples with positive serology with indirect hemagglutination methods, indirect immunofluorescence and other ELISAs, the sensitivity is 99.3%.

b) Specificity: over a panel of 75 samples with negative xenodiagnostic and serology, the specificity is 98.7%. Over another panel of 200 samples with negative serology by indirect hemagglutination methods, indirect immunofluorescence and other ELISAs, the specificity is 100%.

c) Population-based studies: in a population including healthy individuals, donors, chagasic individuals and other pathologies, the correlation with other reference methods was 99.6%.

Despite the sensitivity of the method, its results, as well as those of any other serological method, constitute auxiliary data for the diagnosis. Therefore, the reports should be considered in terms of probability. In this case, more or less probability of *T. cruzi* parasitemia.

Any reactive result must be verified by other technique, such as immunofluorescence, complement fixation, indirect hemagglutination, direct hemagglutination, etc.

Verify any Reactive result using an additional technique. Mind the "Fátala Chaben" Institute recommendations, according to which the immunodiagnosis of the infection should be performed with a minimum of two of the following methods: indirect immunofluorescence, indirect hemagglutination, ELISA, (latex) particles agglutination, validated by the Centro Nacional de Referencia.

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Wiener Laboratorios S.A.I.C.