Chrom Protein C

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BIBLIOGRAFIA

- Walker FJ (1981) Regulation of activated protein C by protein S. The role of phospholipid in factor Va inactivation, J BiolChem, 256: 11128-31.
- Taylor FB, Lockhart MS (1985) A new function for activated protein C: activated protein C prevents inhibition of plasminogen activators by releasate from mononuclear leukocytes--platelet suspensions stimulated by phorbol diester, Thromb Res, 37(1): 155-64.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th ed. CLSI: H21-A5.
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 4th ed., 2001.

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- Coagulation Control N e Coagulation Control P de Wiener lab.

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MÉTODO DE CONTROLE DE QUALIDADE

Coagulation Control N e Coagulation Control P da Wiener lab.

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REFERÊNCIAS

- Walker FJ (1981) Regulation of activated protein C by protein S. The role of phospholipid in factor Va inactivation, J BiolChem, 256: 11128-31.
- Taylor FB, Lockhart MS (1985) A new function for activated protein C: activated protein C prevents inhibition of plasminogen activators by releasate from mononuclear leukocytes--platelet suspensions stimulated by phorbol diester, Thromb Res, 37(1): 155-64.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 4th ed., 2001.

Chrom Protein C

For the quantitative determination of protein C by chromogenic method

SUMMARY

Protein C is a vitamin K dependent protein, together with its cofactor protein S, acts as a coagulation inhibitor, through the regulation of activated factors V and VIII. It is a serine protease which is activated to activated protein C by thrombin-thrombomodulin complex action on the surface of endothelial cells. Deficiency of both inherited and acquired protein C can lead to thrombotic events with varying degrees of severity. The most common manifestation of the disease is deep vein thrombosis of the lower members, complicated or not with pulmonary embolism.

PRINCIPLE

Protein C in plasma is activated by the snake Agkistrodon contortrix venom. The amount of activated protein C is determined by measuring the rate of hydrolysis of a chromogenic substrate, specific for protein C (Pad-Pro-Arg-pNA). The release of p-nitroaniline is measured at 405 nm and is proportional to the level of protein C.

PROVIDED REAGENTS

A. Reagent A: protein C chromogenic substrate (Pad-Pro-Arg-pNA), 6 µmol. Lyophilized.

B. Reagent B: protein C activator containing 0.8 Units/vial snake venom extract. Lyophilized

C. Reagent C: 6.1 g/L buffer Tris, 12.7 g/L NaCl and 0.1% albumin, pH 8.4.

Calibrator 1: protein C standard (see concentration on the label). Lyophilized.

Calibrator 2: protein C standard (see concentration on the label). Lyophilized.

Calibrator 3: protein C standard containing (see concentration on the label). Lyophilized.

NON-PROVIDED REAGENTS

- Distilled water.
- Wiener lab.'s Coagulation Control N and Coagulation Control P.

INSTRUCTIONS FOR USE

Reagent C: ready to use.

Reagent A, B & Calibrators: reconstitute with the distilled water volume indicated on the label. Cover and let stand for 30 minutes. Mix gently without stirring until obtaining a homogeneous solution.

WARNINGS

Reagents are for "in vitro" diagnostic use.

Calibrators were prepared from non-reactive material for HBsAg, antibodies to HCV and antibodies to HIV 1+2. However,

similarly to blood samples calibrators should be handled as if they were infective.

Use the reagents according to the working procedures for clinical laboratories.

All reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents are stable in refrigerator (2-10°C) until the expiration date shown on the box.

After reconstitution they are stable:

Reagent	37°C	RT (< 25°C)	2-10°C	-20°C
Reagent A	8 hours	1 week	1 month	6 months
Reagent B	-	3 days	1 week	6 months
Calibrator 1-3	-	8 hours	2 days	6 months

SAMPLE

Citrated plasma

a) Collection: obtain blood carefully (avoiding stasis or trauma) and place into a tube with anticoagulant in an exact 9+1 proportion (e.g.: 4.5 ml blood + 0.5 ml anticoagulant). Mix gently. Centrifuge for 15 minutes at 2,500 g and separate plasma before 30 minutes.

b) Additives: use Wiener lab's Anticoagulante TP or 130 mmol/l (3.8%) or 109 mmol/l (3.2%) sodium citrate to obtain plasma.

c) Known Interfering Substances: do not use EDTA or heparin to obtain plasma.

No interferences by heparin are observed up to 2 IU/ml.

Refer to Young, D.S. in references for effect of drugs on the present method.

d) Stability and Storage Instructions: plasma must be stored at room temperature until the test is performed. This period should not last for more than 4 hours. In case the sample cannot be processed within this period, the plasma can be frozen up to 6 months at -20°C.

REQUIRED MATERIAL (non-provided)

- Hemolysis tubes
- Pipettes and micropipettes for measuring the stated volumes
- Water bath at 37°C
- Stopwatch
- Spectrophotometer or coagulation analyzer.

ASSAY CONDITIONS

- Wavelength: 405 nm

- Reaction temperature: 37°C
- Reaction time: 8 minutes
- Final Volume: 1.15 mL for endpoint technique and 0.75 mL for kinetic technique.

PROCEDURE

Kinetic technique

1:6 Substrate Buffer: dilute Reagent A 1:6 with Reagent C (1 part Reagent A + 5 parts Reagent C). Preheat at 37°C. In a tube labeled U (Unknown) or C (Calibrator) place:

	U or C
Sample	50 uL
Reagent B	100 uL

Mix and incubate for exactly 5 minutes at 37°C. Then add:

Substrate Buffer 1:6	0,60 mL
-----------------------------	---------

Mix and read the absorbance at 405 nm for 3 minutes ($\Delta A/\text{min}$) recording absorbances every minute. Calculate protein C activity (%) using the absorbance value by interpolation from the calibration curve.

CALIBRATION CURVE

Plot the absorbance obtained for each calibrator ($\Delta A/\text{min}$ or A) according to protein C % for each Calibrator indicated on the label. Protein C concentration in the patient's plasma is determined by interpolation on the calibration curve.

REFERENCE VALUES

70-130 %

Each laboratory should establish its own reference values based on the techniques and devices used.

QUALITY CONTROL METHOD

Wiener lab.'s Coagulation Control N and Coagulation Control P.

PROCEDURE LIMITATIONS

See Known Interfering Substances and Stability and Storage Instructions under SAMPLE.

Recalibration is required for each batch of reagents and for each instrument used.

Heparin values in plasma above 2 IU heparin/mL affect the test results.

PERFORMANCE

a) Reproducibility: reproducibility was determined with various samples (in series and day-to-day). The following results were obtained:

Intra-assay precision

Level	S.D.	C.V.
106.2%	2.23%	2.10%
22.2%	1.32%	5.95%

Inter-assay precision

Level	S.D.	C.V.
106.2%	2.87%	2.70%
22.2%	1.31%	5.90%

b) Measuring range: 10 - 170%

PARAMETERS FOR AUTOANALYZERS

Refer to the specific applications of each autoanalyzer.

WIENER LAB PROVIDES

3 x 1 mL Reagent A
 3 x 1 mL Reagent B
 1 x 60 mL Reagent C
 1 x 1 mL Calibrator 1
 1 x 1 mL Calibrator 2
 1 x 1 mL Calibrator 3
 (Cat. N° 1705013)

REFERENCES

- Walker FJ (1981) Regulation of activated protein C by protein S. The role of phospholipid in factor Va inactivation. J Biol Chem. 256: 11128-31
- Taylor FB. (1985) A new function for activated protein C: activated protein C prevents inhibition of plasminogen activators by releasate from mononuclear leukocytes-platelet suspensions stimulated by phorbol diester. Thromb Res. 37(1): 155-64.
- Clinical and Laboratory Standards Institute (2008) Collection. Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline. 5th edn. CLSI: H21-A5
- Young D.S. - "Effects of Drugs on Clinical Laboratory Tests". AACCC Press. 4th ed.. 2001.

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Wiener Laboratorios S.A.I.C.