Anticoagulantes

Anticoagulante W
Anticoagulante G
Anticoagulante TP

Anticoagulante W

USO PREVISTO

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- 6 frascos conta-gotas x 50 ml (Cód. 1895002).

REFERÊNCIA

- International Committee for Standardization in Hematology (ICSH) - Am. J. Clin. Path. 68/4:505 (1977).

Anticoagulantes

Anticoagulante W
Anticoagulante G
Anticoagulante TP

Anticoagulante W

INTENDED PURPOSE

Reagent for use as an anticoagulant for blood samples in hematology and blood chemistry tests requiring ethylenediaminetetraacetic acid (EDTA), such as red blood cell counts, hemoglobin determination, reticulocyte count, platelet count, hematocrit, erythrocyte sedimentation rate, blood typing and classification, and other chemical determinations requiring plasma (except sodium, potassium, and calcium); for qualitative or quantitative analysis by automated or non-automated methods. Intended for use by biochemical professionals in clinical laboratories, medical clinics, and hospitals, with samples from the general population.

PROVIDED REAGENT

Anticoagulante W: 0.342 mol/l EDTA sodium and potassium salt solutions, pH 7.2.

INSTRUCTIONS FOR USE

Ready to use.

WARNINGS

The Reagent is for "in vitro" diagnostic use.

H319: Causes serious eye irritation. P280: Wear eye protection/ protective gloves/protective clothing. P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P337+P313: If eye irritation persists: Get medical advice/attention. Use the reagent according to the working procedures for clinical laboratories.

Reagent and samples should be discarded according to the local regulations in force.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

STABILITY AND STORAGE INSTRUCTIONS

Anticoagulante W is stable at room temperature (< 25°C) until the expiration date stated on the box.

INSTABILITY OR DETERIORATION OF REAGENT

Its deterioration is indicated by the presence of discoloration or sedimentation.

SAMPLE STABILITY

The blood collected with ethylenediaminetetraacetic acid

(EDTA) shows stability with cellular components and no signs of hemolysis for up to 8 days from collection.

Blood cell counts, hemoglobin, reticulocytes, platelets and hematocrit do not show variations within 24 and 48 hours in blood samples stored at room temperature and at 2-10°C, respectively.

Erythrocyte sedimentation rate and blood for morphologic tests can be performed in blood samples stored at room temperature within 3 to 6 hours from collection or stored for up to 24 hours in refrigerator.

PROCEDURE

One drop (70 µl) inhibits up to 9 ml of blood coagulation. Under such conditions, the dilution error is 0.8%, which is adequate for any routine hematology determination.

For smaller sample quantities, use a reduced amount of anticoagulant. Use 20 µl for volumes up to 2.5 ml and 50 µl to collect up to 7 ml blood.

In case maximum accuracy is required, anticoagulant can be dried at 37°C or at 50°C.

If bone marrow or hematic puncture fluids are used (CSF, ascitic fluids, etc.), use the advised blood proportion.

PROCEDURE LIMITATIONS

The use of higher amounts of anticoagulant may yield falsely decreased hematocrit values.

WIENER LAB. PROVIDES

- 6 dropper bottles x 50 ml (Cat. Nr. 1898552).

Anticoagulante G

INTENDED PURPOSE

Reagent for use as an anticoagulant for blood samples, exclusively for the qualitative or quantitative determination of glucose in the general population, using automated or non-automated methods. Anticoagulant G combines two fundamental principles: an optimal antiglycolytic agent, potassium fluoride, and the anticoagulant of choice, EDTA. Intended for use by biochemical professionals in clinical laboratories, medical clinics, and hospitals.

PROVIDED REAGENT

Anticoagulante G: 0.274 mol/l EDTA sodium and potassium salt solutions, and 0.86 mol/l fluoride, pH 7.2.

INSTRUCTIONS FOR USE

Ready to use.

WARNINGS

The Reagent is for "in vitro" diagnostic use.

H302: Harmful if swallowed. P264: Wash hands and exposed skin thoroughly after handling. P501: Dispose of contents/ container in accordance with local regulation.

Use the reagents according to the working procedures for clinical laboratories.

Reagent and samples should be discarded according to the local regulations in force.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

STABILITY AND STORAGE INSTRUCTIONS

Anticoagulante G is stable at room temperature (< 25°C) until the expiration date stated on the box.

INSTABILITY OR DETERIORATION OF REAGENT

Presence of discoloration or sedimentation indicates deterioration.

SAMPLE STABILITY

Glucolysis is an enzymatic process observed "in vitro" which starts at the moment of sample collection and remains active even in frozen samples.

The range of glucose consumption at 37°C is 0.10 to 0.20 g/l/hour. Leukocytes or bacteria contamination may yield false hypoglycemia in samples stored at room temperature.

PROCEDURE

One drop (70 µl) inhibits up to 9 ml of blood coagulation. Use 20 µl for volumes up to 2.5 ml and 50 µl to collect up to 7 ml blood.

PROCEDURE LIMITATIONS

Note that since fluoride is an enzymatic inhibitor, **Anticoagulante G** should not be used in reactions involving enzymes (e.g. urea with urease).

WIENER LAB. PROVIDES

- 6 dropper bottles x 50 ml (Cat. Nr. 1890552).

Anticoagulante TP

INTENDED PURPOSE

Reagent for use as an anticoagulant in blood samples for the determination of prothrombin time, erythrocyte sedimentation rate, coagulation studies, thromboplastin generation tests, fibrinogen and other coagulation factor determination, platelet count, and other reactions requiring citrated blood, qualitatively or quantitatively, by automated or non-automated methods, depending on the test. Intended for use by biochemical

professionals in clinical analysis laboratories, medical clinics, and hospitals, in samples from the general population.

PROVIDED REAGENT

Anticoagulante TP: 130 mmol/l dihydrated trisodium citrate solution, pH 7.2.

INSTRUCTIONS FOR USE

Ready to use.

WARNINGS

The Reagent is for "in vitro" diagnostic use.

Use the reagent according to the working procedures for clinical laboratories.

Reagent and samples should be discarded according to the local regulations in force.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

STABILITY AND STORAGE INSTRUCTIONS

Anticoagulante TP is stable at room temperature (< 25°C) until the expiration date stated on the box.

INSTABILITY OR DETERIORATION OF REAGENT

Its deterioration is indicated by the presence of discoloration or sedimentation.

SAMPLE STABILITY

Samples collected with citrate concentration of 120 to 150 mmol/l show maximum stability for prothrombin time determination for up to 8 hours after collection. In this condition, blood samples could be stored at 4°C for up to 7 days keeping constant V and VIII coagulation factors' activity.

PROCEDURE

For coagulation tests, the adequate blood/anticoagulant ratio is 9+1. For example, 4 drops of Anticoagulante TP for 2.5 ml blood.

To perform erythrocyte sedimentation rate, the indicated ratio (ICSH) is 4+1. For example, 2 ml blood + 7 drops (0.5 ml) of Anticoagulante TP.

PROCEDURE LIMITATIONS

Note that other anticoagulant/sample or citrate concentration than the one specified will affect Prothrombin Time. Therefore, it is recommended to control the anticoagulant dose when sample is collected.

WIENER LAB. PROVIDES

- 6 dropper bottles x 50 ml (Cat. Nr. 1895002).

REFERENCES

- International Committee for Standardization in Hematology (ICSH) - Am. J. Clin. Path. 68/4:505 (1977).

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Wiener Laboratorios S.A.I.C.